



# PALOMAR LUX1540™ FRACTIONAL LASER

A pictorial description of the skin resurfacing treatment

1



The skin is cleansed of all makeup and laser eye shields are applied.

The skin is retracted as the shields are being applied to ensure the maximal amount of skin is treated.

3



The operator ensures that areas such as the laughter lines are treated thoroughly.

2



Lux1540™ treatments do not require dye application or the use of topical anesthetics. The sapphire tip of the laser is chilled keeping the skin comfortable during the treatment.

The skin is gently stretched as the area is being treated to ensure even coverage. Generally three passes are performed on an area before moving on to the next region.

4



The tip of the Lux1540™ can be applied to the edge of the eye shields allowing for treatment of the skin immediately below the eye and under the eyebrow.



The contact cooling of the Lux1540 ensures that even the delicate skin of the under eye area can be treated.



# PALOMAR LUX1540™ FRACTIONAL LASER

5



Another popular area for Lux1540 treatments are the lip lines and the area around the mouth.

6



Fractional non-ablative skin resurfacing can be applied to any area of the body, including the face, neck, chest, arms and hands.

7



The skin after a full face Lux1540™ treatment showing slight pinkness and faint swelling. The skin is not broken and an application of light camouflage will completely conceal any residual effects.

## WRINKLES AROUND THE EYE

BEFORE



Upper and lower eye areas and laughter lines.

AFTER



Visible smoothing and reduction of wrinkles with reduced dark circles.

## SURGICAL SCAR

BEFORE



Surgical scar resulting from excision of birthmark in childhood.

AFTER



Five weeks after third Lux1540™ treatment.